



Yuzhpolymetal-Holding Group

# Guard-SRS

— stationary Raman analyzer for rapid identification of bottled liquids



- \* Dimensions: 500×330×1400 mm
- \* Allowable objects of analysis:
  - glass or plastic bottles**
  - or other containers**
  - in volume from 100 ml to 5 liters**
- \* Analysis time: 1 – 10 s
- \* Identified chemicals:
  - ✓ **liquid explosives and their components**
  - ✓ **caustic and corrosive liquids**
  - ✓ **flammable liquids**

## Main advantages:

- Wide range of identified chemicals
- Long service life
- No expendable materials
- The method of analysis, safe for humans and the environment

The **Guard-SRS** analyzer is optimized for use in the preflight inspection of carry-on baggage and luggage of passengers, allowing to identify all the basic substances of the list of dangerous substances prohibited for carriage on board (ICAO Doc 9284-AN/905)

“Yuzhpolymetal-Holding” Group

✉ 56, Varshavskoye street, Moscow 117638 Russia

☎ +7 499 613-11-77, +7 499 317-31-55, +7 499 317-31-66

e-mail: [analizator@list.ru](mailto:analizator@list.ru)

[www.analizator.ru](http://www.analizator.ru)